DEVICE DESCRIPTION

This non-sterile volar plate system offers several sizes of volar T plates made from stainless steel (ASTM F-138). The distal head of the plate contains two rows of fixed angle screw holes that accept 2,7 mm screws. Radial and ulnar sheaths are secured into the plate recesses to cover the screw heads and to lock the screws in the fixed orientation. The 2,7mm screws vary in length from 10mm to 34mm. The proximal end of the plate is secured to the bone with 3,3mm screws. The 3,3mm screws vary in length from 8mm to 18mm. K-wires may also be used for stabilization of bone fragments.

INDICATIONS FOR USE

This device is indicated for temporary or permanent volar fixation and stabilization of fractures and osteotomies involving the distal radius.

This is single use device is for cementless use only.

CONTRAINDICATIONS

This device is contraindicated in the following situations:

- Lack of sufficient sound bone to seat the screws including that due to skeletal immaturity or osteoporosis.
- Metal allergies or sensitivity.
- Infection at or near the site of implantation.
- Distant infection.
- Fractures exceeding the length of the device.

WARNINGS AND PRECAUTIONS

The Coverloc™ Volar Plate System implants are for single use only and should never be resterilized after contact with body tissues or fluids in order to prevent any risks of cross-contamination. Do not implant plates and screws of dissimilar materials. Do not implant plates and K-wires of dissimilar materials.

Implanting internal fixation devices of dissimilar metals may result in electrolytic corrosion at the interface of the devices.

Do not overbend or rebend the plate. Overbending or rebending the plate will reduce its mechanical strength and may result in device breakage.

The Coverloc™ Volar Plate System implants have not been evaluated for safety and compatibility or tested for heating or migration in the MR environment.

Although published studies indicate that the MRI procedure has minimal effects on the material in the Coverloc™ Volar Plate System implants, patients should be informed that there are several different manufacturers and generations of MRI equipment available, and Tornier Inc. cannot make claims regarding the compatibility of Coverloc™ Volar Plate System implants with any specific MRI unit. It is recommended that patients and physicians consult with the MRI equipment manufacturers to discuss the compatibility of the Coverloc™ Volar Plate System implants with the MRI equipment before undergoing any test.

For any questions, contact . Alternatively you can directly contact the manufacturer.

CLEANING AND DISINFECTING

WARNING: The following Cleaning and Disinfecting instructions have been validated. Failure to follow all steps may result in an improperly cleaned and sterilized device (Non-Sterile). See **HOW SUPPLIED** for sterilization instructions.

- PREPARATION FOR CLEANING**
 - Remove dried on soil from the instrument, particularly in areas such as mechanical joints or crevices.
- MANUAL CLEANING**
 - Clean and remove all foreign debris prior to disinfecting and sterilizing.
 - CAUTION: Tornier has not validated the use of an automated washer-disinfector for cleaning, disinfecting or sterilizing the device.
 - Immediate rinsing and cleaning after use with an enzymatic or alkaline cleaning detergent will effectively remove and prevent drying of adherent blood, mucus, etc. Prepare solutions per manufacturer instructions. Cleaning solutions include, but are not limited to: ENZOL™ enzymatic, Neodisher® Mediclean Forte, Alconox®Terg-A-Zyme™ and Thermosept® Alka Clean.
 - Scrub instrument with a soft brush, paying particular attention to areas where debris might accumulate. Always avoid any harsh materials that can scratch or mar the surface of the instrument.
 - Rinse the instrument thoroughly with distilled water or preferably demineralized sterile water.
 - Inspect instrument for visible soil. Repeat cleaning if soil is visible and re-inspect.
- ULTRASONIC CLEANING**
 - Place instrument in an ultrasonic cleaner for 20 minutes and process per the ultrasonic cleaner unit's instructions in a pH neutral solution.
 - Rinse the instrument thoroughly with water following the ultrasonic process.
 - Inspect instrument for visible soil. Repeat cleaning if soil is visible and re-inspect.
- MANUAL DISINFECTION**
 - Instruments MUST be cleaned prior to disinfection, as blood albumin will impair the bactericide effectiveness of the solution.
 - CAUTION: Tornier has not validated the use of an automated washer-disinfector for cleaning, disinfecting or sterilizing the device.
 - Immerse the instrument in disinfection solutions for a minimum of 20 minutes. Suitable disinfection solutions can include, but are not limited to: CIDEX™, WAVICIDE™-01, Gigasept®, Kohrsolin®, and equivalent solutions. Use the supplier's instructions for preparing the solution. CAUTION: mercurial solutions are not recommended, as they corrode metal parts.
 - Manipulate all moving parts a minimum of 5 times while completely immersed.
 - After disinfection, rinse the instrument with distilled water or preferably demineralized sterile water. Repeat rinse if any residues are visible.
- AUTOMATIC WASHING**
 - CAUTION: Tornier has not validated the use of an automated washer-disinfector for cleaning, disinfecting or sterilizing the device.
 - A cleaning process done out of qualification ranges can lead to sterility or toxicity issues.

POSSIBLE ADVERSE EVENTS

As with all implantable devices, risks can include bone resorption, poor bone formation, tissue or nerve damage as the result of surgical trauma, allergic reaction to the materials used in the device, infection, screw pullout, device fracture, device loosening, and/or fracture nonunion. Adverse events are more likely if the device is not removed after fracture healing occurs.

DIRECTIONS FOR USE WITH THE VOLAR SURGICAL APPROACH USING THE COVERLOC™ VOLAR PLATE: (excerpts taken from Volar Plate with CoverLoc™ Technology Surgical Technique, #30009)

- Make a 3-6cm longitudinal incision along radial border of the FCR tendon beginning at the wrist flexion crease and proceeding proximally.
- The pronator quadratus is exposed and released sharply from its radial insertion in a region which maintains good tissue on either side for later repair.
- Additionally, the insertion of the brachioradialis can be released, eliminating an important deforming force on the distal fragment.
- Excise all interposed soft tissue and perform a preliminary reduction of the fracture.
- Select the appropriate plate size and attach the corresponding drill guide on the back table.
- Determine the correct position for the plate by placing the plate along the volar surface of the radius, and aligning it distally along the distal radius margin and proximally along the radial shaft. Evaluate fluoroscopically.
- Note: The use of provisional fixation is provided for through k-wire holes in the proximal stem of the plate.
- Use the 2,5mm drill and guide to drill a hole through the center of the oval slot.
- Measure carefully using the depth gauge for the cortical (proximal) screw. Insert the 3,3mm cortical screws using the hex driver. This cortical screw placement allows proximal and distal plate adjustment as may be indicated.
- With the fracture held reduced, insert the 2,0mm drill sleeve into the desired hole until it is fully seated. The #2 or #3 hole are generally preferred. NOTE: The sleeve does not screw into the hole.
- Evaluate screw angle under fluoro before drilling. If desired, K-Wire drill sleeve may be used first and a 0,045" K-wire inserted to evaluate screw angle.
- Drill through both cortices using the drill sleeve with the 2,0mm drill bit.
- Measure length carefully, using the metaphysis (distal) depth gauge. NOTE: The CoverLoc™ technology allows the screw to continue turning in the hole to lag the fragment to the plate (there are no threads on the screw head). Due to the lagging capability, ALWAYS select the shorter length if between sizes.
- Select the proper length screw from the caddy. The #8 star drive will "pick and stick" from the caddy.
- Insert the 2,7mm screw until the line on the driver shaft reaches the recess on the drill guide. This line provides a preliminary seating and final tightening is done after screws are in place and guide is removed.
- Repeat the process of inserting the 2,0 mm drill sleeve into a hole: drilling, measuring, and inserting the distal screws.
- Upon completion of fixation, fluoroscopic evaluation is essential to determine fracture reduction and radius alignment. Care must also be taken to ensure that no screw protrudes beyond the dorsal cortex to avoid tendon irritation or rupture. NOTE: Be cognizant of the prominence of Lister's Tubercle dorsally; it may hide a prominent screw in the 3rd extensor compartment. True lateral and lunate fossa views of the distal radius are also evaluated to ensure that no screw has penetrated the distal radius articular surface.
- Drill, measure, and insert remaining proximal screws that may be required for adequate fixation.
- Ensure all screw heads are fully seated before placing the CoverLoc™ sheath covers.
- Fix CoverLoc™ screws using a firm finger-tip grip. If the CoverLoc™ does not seat fully, remove and fully seat each screw.
- The pronator quadratus is reapproximated over the plate with interrupted non-absorbable sutures.
- Evaluate final construct under fluoroscopy.
- Close the surgical incision in the normal fashion.

HOW SUPPLIED

This device is supplied **NON-STERILE** and must be **CLEANED** (see CLEANING AND DISINFECTING) and **STERILIZED** prior to use. Steam sterilization at 134°C (273°F) for 30 minutes is recommended. Flash sterilization is not recommended.

STORAGE INSTRUCTION

Store in a cool dry place and keep away from direct sunlight.

	CATALOG NUMBER		DO NOT USE IF PACKAGE IS DAMAGED
	BATCH CODE		NON-STERILE
	AUTHORIZED REPRESENTATIVE OF THE EUROPEAN COMMUNITY		CAUTION: FEDERAL LAW (USA) RESTRICTS THIS DEVICE TO SALE BY OR ORDER OF A PHYSICIAN
	CAUTION		CE MARK & IDENTIFICATION NUMBER OF NOTIFIED BODY
	CONSULT INSTRUCTIONS FOR USE		MANUFACTURER
	NOT MADE WITH NATURAL RUBBER LATEX		

	Tornier, Inc. 10801 Nesbitt Avenue South Bloomington, MN 55437 USA Tel: (952) 921-7100 Fax: (952) 236-4007		Tornier SAS 161 Rue Lavoisier 38330 Montbonnot Saint Martin France Tel : +33(0) 476 613 500 Fax : +33(0) 476 613 533
---	---	---	---

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