Stabilis™ Ankle Arthrodesis Plate System

19-5113 Rev D

DESCRIPTION

The Stabilis™ Ankle Arthrodesis Plate System consists of minimally contoured plates, cover sheaths, screws and washers made of stainless steel, with instruments to facilitate fixation in ankle arthrodesis and fracture fixation, fusion, and osteotomy of the bones of the foot.

The plates are not differentiated for right or left sides and include an Anterior Tibio-Talar Plate (two sizes), a Lateral Tibio-Talar Plate (two sizes), a Lateral Tibio-Talar-Calcaneal Plate and an optional augmentation plate. The screws are provided in one diameter and in various lengths. The washers are provided in one size.

INDICATIONS FOR USE

The Stabilis™ Ankle Arthrodesis Plate System is intended for arthrodesis of the ankle joint and distal tibia, fractures, osteotomies, and fusions of small bones including the foot and ankle.

CONTRAINDICATIONS

The Stabilis™ Ankle Arthrodesis Plate System is contraindicated for the following conditions:

- Sepsis, systemic infection, elevated WBC count, fever and/or local inflammation.
- Complete talar necrosis or inadequate bone stock to permit stabilization of the arthrodesis.
- Persisting skin lesion or poor skin coverage around the ankle joint that would make the procedure unjustifiable.
- Neuromuscular or psychiatric disorders which might jeopardize fixation and post-operative care.
- Suspected or documented metal allergy or intolerance.
- Inadequate neurovascular status.
- Severe longitudinal deformity.

WARNINGS (See also the Patient Counseling Information Section)

Strenuous loading, excessive mobility, and articular instability all may lead to accelerated wear and eventual failure by loosening, fracture, or dislocation of the device. If excessive loading cannot be prevented, an implant should not be used. Abnormal or excessive forces could lead to delayed union, non-union, or failure of the implant. Patients should be made aware of the increased potential for device failure if excessive demands are made upon it.

Implants are mechanical devices that can be worn away, fatigued or broken. An implant site may become infected, painful, swollen, or inflamed. The status of the adjacent bone and soft tissue may be inadequate to support the implant, or may deteriorate in time resulting in instability, deformity, or both. The benefits from implant surgery may not meet the patient’s expectations or may deteriorate with time, necessitating revision surgery to replace the implant or to carry out alternative procedures. Revision surgeries with implants are not uncommon.

The Stabilis™ Ankle Arthrodesis Plating System has not been evaluated for safety and compatibility in the MR environment. The system has not been tested for heating or migration in the MR environment.

PRECAUTIONS

Implants should only be handled with blunt instruments to avoid scratching, cutting or nicking the device.

Meticulous preparation of the implant site and selection of the proper size implant increase the potential for successful outcome.

This implantable product is for single use only. An implant should never be re-sterilized after contact with body tissues or fluids.

ADVERSE EVENTS

Potential adverse events reported with bone plates and screws include pain, loosening, fracture, non-union, misalignment, dislocation, or infection or swelling at the implant site. There have been some reports of patients with metal sensitivity reactions following plate and screw implantation. Implantation of implant materials may result in foreign body reaction adjacent to the implant site. Injury to the surrounding bone, nerves, blood vessels, tendons, or soft tissues can occur as a consequence of implanting this device.

SURGICAL PROCEDURES

A surgical procedure for use of this implant device is available. It is the responsibility of the surgeon to be familiar with the procedure before use of these products. In addition, it is the responsibility of the surgeon to be familiar with relevant publications and consult with experienced associates regarding the procedure before use.

PATIENT COUNSELING INFORMATION (See also Warnings)

Providing each patient scheduled for implant surgery with documented counseling of potential complications and alternatives, which may include non-implant procedures such as soft tissue reconstruction, prior to surgery is necessary. In addition to the patient related information contained in the Warnings and Adverse Events sections, the following information should be conveyed to the patient:

- Adverse effects may necessitate reoperation, revision, or fusion of the involved joint.
- While the expected life of an implant is difficult to estimate, it is finite. These components are made of foreign materials which are placed within the body. However, due to the many biological, mechanical and physiochemical factors which affect these devices, the components cannot be expected to withstand the activity level and loads of normal healthy bone for an unlimited period of time.

CLEANING AND STERILIZATION INFORMATION WARNINGS

- Clean, not sterile, implants and instruments must be sterilized prior to each use.
- The most common cause of instrument breakage is misuse. Specialty instrumentation should never be used for tasks it was not specifically designed to perform. Misuse of an instrument may result not only in damage to the instrument, but also trauma to the patient and operating room personnel.
- If not handled properly, instruments with cutting edges or sharp corners may compromise sterility by tearing surgical gloves. If an instrument tip becomes bent, chipped, or otherwise damaged the instrument should be replaced or repaired before further use. Attempts to straighten bends are not advised as the metallurgical integrity of the metal may be compromised in the process, and the instrument may subsequently break during use.

CLEANING INSTRUCTIONS

New implants and instruments are packaged clean, not sterile and should be assumed to be contaminated. Implants and instruments must be cleaned and sterilized before each use.

For your safety, be familiar with the procedures for handling contaminated materials at your facility before following these instructions.

Clean instruments as soon as possible after use and avoid allowing soiled instruments to dry. Immerse, use damp towels or sponges with deionized or distilled water to keep soiled instruments moist prior to cleaning. Manually scrub with a clean, soft-bristled brush or soft cloth in mild detergent following the detergent manufacturer’s instructions for use. Avoid using extreme detergent concentration levels. Detergents with enzyme cleaners or separate enzyme solutions may be used, but their effectiveness has not been evaluated. Warm or hot water should be used to aid in cleaning. Ultrasonic cleaning is recommended for instruments with internal surfaces or crevices which may be hard to clean manually. pH neutral cleaners (pH 6.0-8.0) are recommended for longer life of the instruments. If acidic or alkaline solutions are used follow the manufacturer’s recommendations for neutralizing the pH by rinsing with water or neutralizing agent. Highly alkaline or acidic cleaners (used in some mechanical washers) are not recommended as they will reduce the life of the instruments and may affect instrument performance. Avoid prolonged exposure to acidic or alkaline solutions and solutions containing chlorides, bromides or iodine.

Use of water soluble lubricants is recommended for instruments with moving parts or those intended to be assembled intraoperatively to another instrument. Rinse thoroughly with deionized water or distilled water. Dry completely before sterilizing. Inspect each instrument and implant thoroughly for cleanliness, especially in broach teeth, internal surfaces and crevices. Instruments should be re-washed if there are any signs of foreign matter or residue. Check instruments thoroughly for damage. Do not use damaged instruments as they may compromise the surgical outcome. Replace damaged instruments before next use.

Use of mechanical washers has not been evaluated by the manufacturer. Qualification of specific wash cycles and equipment should be completed by the user.

STERILIZATION


The recommended sterilization method is steam sterilization. The following sterilization cycles have been shown to produce a sterility assurance level of 10⁻⁶ when parts have been cleaned to the instructions above. Other similar steam cycles and cleaning procedures may be used but have not been evaluated. Sterilization qualification was performed using specific equipment and procedures. Use of cycles, equipment and procedures other than those listed should be qualified by the user. Do not exceed 300°F/149°C.

Wrapped Pre-vacuum Steam Sterilization 4 minutes minimum @ 270-275°F (133-135°C). Drying time: 30 minutes minimum.

After sterilization, remove the materials from their packaging or the sterilization tray using accepted sterile technique. Ensure that implants are at room temperature prior to implantation.

LIMITED WARRANTY

Tornier, Inc. warrants that this product meets the manufacturer’s specifications and is free from manufacturing defects at the time of delivery. This warranty specifically excludes defects resulting from misuse, abuse, or improper handling of the product subsequent to receipt by the purchaser. Tornier, Inc. does not warrant the outcome of the surgical procedure.

	CATALOG NUMBER
	BATCH CODE
	CAUTION
	CONSULT INSTRUCTIONS FOR USE
	DO NOT REUSE
	NOT MADE WITH NATURAL RUBBER LATEX
	DO NOT USE IF PACKAGE IS DAMAGED
	NON-STERILE
	CAUTION: FEDERAL LAW (USA) RESTRICTS THIS DEVICE TO SALE BY OR ORDER OF A PHYSICIAN
	MANUFACTURER

	Tornier, Inc. 10801 Nesbitt Avenue South Bloomington, MN 55437 USA Tel: (952) 921-7100 Fax: (952) 236-4007
-------------------------------------------------------------------------------------	---------------------------------------------------------------------------------------------------------------------------

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