**TORNIER** 

**EN**

**INSTRUCTIONS FOR USE**  
**AFFINITI™ SHOULDER PROSTHESIS**  
BAW-0759 Rev F

**IMPORTANT:** The implantation of a joint prosthesis requires knowledge of anatomy, biomechanics, and reconstructive surgery of the musculoskeletal system and may be performed only by a qualified surgeon. The surgeon must operate in accordance with current information on the state of scientific progress and the art of surgery. Patients should be aware of the possible complications that can occur as a result of disregarding the precautions listed below.

**Caution:** Federal (USA) law restricts this device to sale by or on the order of a physician.

**DESCRIPTION**  
The Tornier Affiniti™ Shoulder System is a non-constrained glenohumeral prosthesis intended for use as a total or hemi-shoulder replacement system. The total shoulder consists of a metal humeral stem, a metal humeral head and an ultrahigh molecular weight polyethylene glenoid. The hemi-shoulder consists of a metal humeral stem and metal humeral heads. The metal humeral heads are offered in following styles: standard, eccentric or extended humeral heads. The humeral stem is available with and without porous coating. For a more detailed description of the implants and their utilization, please refer to technical documentation, or contact your Tornier representative. It is essential to implant the Affiniti™ Shoulder System with the Affiniti™ instrumentation specifically designed for this purpose. Affiniti™ implants are not compatible with components from other manufacturers. Glenoid components are labeled "for cemented use only" and are indicated only for use with bone cement. Humeral stems are indicated for press-fit uncemented use or for use with bone cement.

**Note:**The humeral stems with porous coating and extended humeral heads (EH)are cleared for distribution in the U.S. only.

#### MATERIALS

The material used in the manufacture of the Affiniti™ implants are as follows: Humeral stems with porous coating are made of titanium alloy (Ti-6Al-4V) according to ASTM F-136 and a sintered titanium bead coating according to ASTM F-1580. Humeral stems without porous coating are made of titanium alloy (Ti-6Al-4V) according to ASTM F-136. The humeral head is manufactured from Cobalt Chrome alloy per ASTM F-1537 (equivalent to ISO 5832-12, "Implants for surgery – Metallic materials – Part 12: Wrought cobalt-chromium-molybdenum alloy") Or ASTM F1058 (equivalent to ISO 5832-7, "Implants for surgery – Metallic materials – Part 7: Forgeable and cold-formed cobalt-chromium-nickel-molybdenum-iron alloy). Glenoid components are made of ultra-high molecular weight polyethylene (UHMWPE) according to ASTM F-648 or ISO 5834-2.

#### INDICATIONS FOR USE

The Affiniti™ Total and Hemi-Shoulder System are indicated for:

- A severely painful and/or disabled joint resulting from osteoarthritis, traumatic arthritis or rheumatoid arthritis;
- Fracture/dislocations of the proximal humerus; where the articular surface is severely comminuted, separated from its blood supply or where the surgeon's experience indicates that alternative methods of treatment are unsatisfactory;
- Other difficult clinical problems where shoulder arthrodesis or resection arthroplasty are not acceptable (e.g., revision of a failed primary component).

The Affiniti™ Hemi-Shoulder System is also indicated for:

- Ununited humeral head fractures;
- Avascular necrosis of the humeral head;
- Rotator cuff tear arthropathy.

**Note:** Humeral stems and glenoids labeled "for cemented use only" are indicated only for use with bone cement. Humeral stems are also indicated for press-fit uncemented use or for use with bone cement.

**Note:**This is a single use device.

**Note:**The humeral stems with porous coating and extended humeral heads (EH) are cleared for distribution in the U.S. only.

#### CONTRAINDICATIONS

**For Total Shoulder or Hemi-Shoulder:**

This device is contraindicated in the following situations:

- Lack of sufficient sound bone to seat and support the implant, including that due to skeletal immaturity or osteoporosis.
- Metal allergies or sensitivity.
- Infection at or near the site of implantation.
- Distant or systemic infection.

**For Total Shoulder Only:**

- Lack of sufficient sound muscle or rotator cuff to seat and support the implant.

#### ADVERSE EFFECTS

As is common with all implantable devices, risks can include component loosening, dislocation, subluxation, bone resorption, poor bone formation, tissue or nerve damage as a result of surgical trauma, infection and/or those events common to any surgery (pulmonary embolism, stroke, heart attack, hematoma, pneumonia, slow wound healing). Additionally, iatrogenic fracture and traumatic fracture below the humeral component have been reported following use of similar devices. Metal sensitivity in patients following total joint replacement have been rarely reported. In the case of revision, all events are more likely and more severe, including difficulty in removing the prior implant and bone cement, challenges in placing the incision and hematoma, as well as wound healing.

#### WARNINGS AND PRECAUTIONS

- Do not re-use this device.
- Do not use with devices from other manufacturers.
- Do not use trials with implants.
- Do not alter or modify the implant.
- Due to the potential interference between the collar of the stem and the skirt on the underside of the head a small number of *eccentric* head and stem combinations may not be possible. In this case, a centered head of the same size may be used.
- If cement is utilized a distal cement restrictor should be placed 1–2cm beyond the tip of the humeral stem to prevent cement extravasations into the distal humeral canal.

#### CAUTIONS

The functional life expectancy of this device has not been determined.

If the surgeon determines this implant is appropriate for a patient with any of the following conditions, the surgeon must advise the patient of the strength limitations of the implant, the possibility of implant dislocation, migration, loosening or failure and the precautions the patient must take to minimize the possibility of an adverse clinical result.

These conditions may impose severe loading on the device:

- Obesity
- Engaging in manual labor or activities
- Certain sports activities
- Patient youth and activity level
- Medical conditions or balance impairments that could lead to falls
- Alcohol or drug addiction
- Mental attitudes or disorders that could result in failure to follow surgeon orders

These conditions may adversely affect the fixation of shoulder implants, if the conditions create a poor environment or poor bone stock:

- Marked osteoporosis, tumors or generally poor bone stock
- Metabolic or systemic disorders, pharmacological treatments
- History of general or local infections
- Severe deformities
- Allergic reactions to implant materials
- Tissue reactions to implant corrosion or wear debris
- Disabilities of other joints

Using a glenoid prosthesis in patients with cuff tear arthropathy could increase the risk of glenoid loosening due to the proximal migration and non anatomic loading.

The Affiniti™ implants have not been evaluated for safety and compatibility or tested for heating or migration in the MR environment.

Although studies indicate that the MRI procedure has minimal effects on material in the Affiniti™ System implants; Patients should be informed that there are several different manufacturers and generations of MRI equipment available, and Tornier Inc. cannot make claims regarding the compatibility of Affiniti™ System implants with any specific MRI unit. It is recommended that patients contact the surgeon or the manufacturer of the MRI equipment to discuss the compatibility of the Affiniti™ System implants with the MRI equipment before undergoing any test.

#### SURGICAL PROCESS

- Pre-operatively:*

The surgeon must be fully conversant with all the aspects of the surgical technique and know the indications and contra-indications of this type of implant. The surgeon must have acquainted himself before the operation with the specific operative technique of the product which is available from the manufacturer. As part of the pre-operative examination, the surgeon must check that no biological, biomechanical or other factors are present that will affect the correct conduct of the operation and the postoperative period. An appropriate range of sizes must be available at the time of the operation. Check that the implant package has not suffered from any deterioration.

- Per-operatively:*

The correct selection of the type and size of the implant appropriate to the patient and the position of the implant are extremely important. Handle the implant using proper aseptic technique including the use of sterile surgical gloves. Leave the protective cover on until the device is needed. The functional surfaces of the implant must not suffer any damage, abrasion or other deterioration. Do not increase the anteverision of the glenoid component as it may result in instability and dislocation. Before closing the incision, clean out all ectopic bone, extraneous cement and bone chips as debris left in the site may cause dislocation, pain, restricted movement or pre-mature wear.

- Post-operatively:*

Patients should be informed about the precautions they must take in everyday life to help maximize the service life of the implant. It is recommended that a regular postoperative follow-up be undertaken to detect early signs of wear, loosening of the prosthesis, etc., and to consider the action to be taken. Normal wear of the implant according to the state of knowledge at the time of its design cannot in any way be considered to constitute a dysfunction or deterioration in the characteristics of the implant. A suitable rehabilitation program must be designed and implemented.

#### STORAGE AND HANDLING

The prosthetic components must be handled and stored with care in accordance with the provisions of the ISO standard 8828. The implant must be stored in their sealed packaging of origin.

#### PACKAGING AND STERILIZATION


The implants are individually packaged and supplied sterile. All implants have been sterilized by gamma irradiation. The expiration date for sterilization must be checked. Only those products implanted before the end of the expiration date may be considered sterile. The packaging and labelling must be checked for integrity. Reject any implant if the packaging is damaged. Every precaution must be taken to ensure sterility when opening the packaging of the implant and when inserting it. Never re-sterilize and implant.

#### RETRIEVAL AND ANALYSIS

Retrieval and Analysis of prosthetic components should be performed in accordance with the provisions of the ISO standard 12891.

#### LIMITED WARRANTY

Tornier, Inc. warrants that this product meets the manufacturer’s specifications and is free from manufacturing defects at the time of delivery. This warranty specifically excludes defects resulting from misuse, abuse or improper handling of the product subsequent to receipt by the purchaser. Tornier, Inc. does not warrant the outcome of the surgical procedure.

<b>Symbols Key:</b>	
	Do not reuse
	Do not use if package is damaged
	Use by
	Sterilization by radiation
	Batch code
	Catalog number
	Caution, consult accompanying documents
	Consult instructions for use
	Manufacturer

	<b>Manufacturer:</b> Tornier, Inc. 10801 Nesbitt Avenue South Bloomington, MN 55437 USA Tel: (952) 921-7100 Fax: (952) 236-4007
---	---

© Copyright 2009 Tornier, Inc.  
Affiniti™ is a trademark of Tornier, Inc.

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Affiniti™ es una marca de Tornier, Inc.