TORNIER 

EN

SOLOSTITCH SYSTEM

CAW-2747 Rev. E

Please read this information carefully before using the Tornier® Solostitch System.

DESCRIPTION
The Tornier® Solostitch System is designed to facilitate arthroscopic repair of a torn rotator cuff with soft tissue augmentation products such as Conexa™ graft.

INDICATIONS
The Tornier® Solostitch System is intended for use in arthroscopic procedures involving soft tissue repairs, where placement of a surgical mesh or patch is desired.

CONTRAINDICATIONS
None

WARNINGS
None

- PRECAUTIONS**
- Do not use this product without reading and understanding these instructions.
 - Do not resterilize. Discard all open and unused instruments.
 - The product is sterile if the package is unopened and undamaged. Do not use if the package seal is broken.
 - Discard the product if mishandling has caused possible damage or contamination.
 - Aseptic technique must be adhered to throughout the procedure.
 - Discard instrument kit if past Use-by-date of the product [as indicated as year (4 digits) and month (2 digits)]. The product is not to be used after the last day of the month indicated.

DIRECTIONS FOR USE
MEASURING PROBE

- While holding the end of the probe without marks, insert the opposite end into the joint cavity.
- Marks on the probe are spaced according to the spacing of the Solostitch™ (½ cm). Every third mark is bold to denote 1 cm increments.
- Measure the anatomical landmarks by counting probe marks and multiplying by 3 ½ mm.
- The probe may be slightly bent to facilitate navigating tissue or for measuring curved dimensions. Care must be used to ensure the device is not excessively bent, substantially weakening the device.

SOLOSTITCH™ TEMPLATE / FREE NEEDLE

- Place a graft that has been cut to the appropriate size and shape in the Solostitch template.
- Snap the template closed.
- Sutures that have been passed through anatomical structures and exit through a cannula can be passed using the suture needle provided in the kit.
- Slots in the Solostitch template are spaced at ½ cm. Cut-outs along the length of the slots are also spaced at ½ cm. This spacing corresponds to the laser marking on the measuring probe.
- Using the slot spacing and the cut-outs along the length of the slot as guides, pass sutures as required by patient anatomy and the desired reconstruction pattern.
- Once the entire suture passing has been completed, unlatch the template and remove the graft from the instrument. Take care to ensure that no sutures are pulled through graft while removing it from the template.

STERILITY
Tornier® Solostitch System is sterilized. Product is sterile unless package has been opened or damaged. **DO NOT RESTERILIZE.** Discard opened or unused containers.

STORAGE AND HANDLING
Store in a cool, dry place and in a manner that protects the integrity of the package and sterile barrier.

PACKAGING AND LABELING
1. Tornier® Solostitch System should be accepted only if the factory packaging and labeling arrive intact.
2. Contact Customer Service if package has been opened or altered.

CAUTION: Federal law (USA) restricts this device to sale by or on the order of a physician.

SYMBOLS KEY

	Do not reuse
	Do not use if package is damaged
	Use by
	Sterile
	Batch code
	Catalog number
	Caution, consult accompanying documents
	Consult instructions for use
	Manufacturer
	Authorized representative in the European Community
	Product contains no latex

CUSTOMER SUPPORT
For product complaints and potential adverse events, please contact your local Sales Representative or Customer Support.

 **Tornier, Inc.**
10801 Nesbitt Avenue South
Bloomington, MN 55437
USA
Tel: (952) 921-7100
Fax: (952) 236-4007

 **Tornier**
161 Rue Lavoisier
38330 Montbonnot Saint Martin
FRANCE
Tel: 33 (0) 476 613 500
Fax: 33 (0) 476 613 533

Solostitch™ and Conexa™ are trademarks of Tornier, Inc.

Tornier® is a registered trademark of Tornier, Inc. These products are covered by one or more U.S. and international patents pending. All rights reserved. Printed in the USA.

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© Copyright 2011 Tornier, Inc.
Tornier, Inc. 10801 Nesbitt Avenue South Bloomington, MN 55437 USA Tel: (952) 921-7100 Fax: (952) 236-4007
Solostitch™,e Conexa™ zijn handelsmerken van Tornier, Inc.