## INSTRUCTIONS FOR USE

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CAW-2847 Rev. D

**IMPORTANT:** The implantation of a joint prosthesis requires knowledge of anatomy, biomechanics and reconstructive surgery of the musculo-skeletal system and may be performed only by a qualified surgeon. The surgeon must operate in accordance with current information on the state of scientific progress and the art of surgery. The patient must be informed of everything related to him and contained in the present instructions for use.

**Caution:** The Federal (United States) Law restricts this device to sale, distribution and use by or on the order of a physician.

**1. Description:** The Latitude Total Elbow prosthesis is provided with two systems: The Latitude and Latitude EV. Both systems are modular and consist of a humeral spool, humeral stem, ulnar stem, ulnar cap, radial stem and head. The systems share the same humeral spools, ulnar cap, and radial components. The Latitude Elbow System is comprised of **un-coated** humeral stems and **un-coated** ulnar stems. The Latitude EV Elbow System is comprised of **coated** humeral stems and **coated** ulnar stems.

The Latitude and Latitude EV systems are intended for use as a total elbow replacement that can be either non-constrained or semi-constrained. The Latitude system is also intended for use as a hemi elbow replacement with anatomic humeral spools.\*

The humeral implant consists of three components: humeral spool, humeral stem, and humeral screw.

The humeral spool in the Latitude system and Latitude EV system is available in centered, offset, anatomic versions.\*

**Note: Both Latitude and Latitude EV systems share the same humeral spools but not the same humeral screws.**

The humeral stems associated with a humeral spool are articulated with the ulnar and radial components.

The humeral stems associated with an anatomic humeral spool are articulated with native ulnar and radial bones.\*

**\*Note: Anatomic and offset humeral spools are not available for sale within the United States.**

The humeral screw is packaged with the humeral stem and assembles the humeral spool and anatomic spool to the humeral stem.

Sterile humeral screws are only available for sale with the Latitude system and outside of the United States.

The ulnar implant consists of two types: non-constrained and semi-constrained ulnar stems.

The ulnar stem without the cap allows for a non-constrained prosthesis.

The ulnar stem with the cap installed creates a semi-constrained prosthesis.

For a more detailed description of the implants, refer to the technical documentation, surgical technique, or contact your Tornier representative. It is essential to implant the Tornier Latitude & Latitude EV Prosthesis with the Tornier Latitude & Latitude EV instrumentation specifically designed for this purpose. Tornier implants must be assembled using Tornier components defined as being compatible with one another.

All the components are not cleared for use in all countries. Please contact your Tornier representative for information about their availability.\*

**Notes:**

The Latitude EV humeral stems and ulnar stems have the word **“COATED”** on the product label description.

• Symbols can be used to identify some implants (labeling or marking). They have the following meaning:

**Left** = Left; **Right** = Right; **SM** = Small; **SM+** = Small+; **ME** = Medium; **ME+** = Medium+; **LA** = Large; **LA+** = Large+.

All the components are not cleared for use in all countries. Please contact your Tornier representative for information about their availability.

**2. Latitude and Latitude EV Implant Compatibility:** Implant compatibility for Latitude and Latitude EV are explained in the respective surgical techniques UCLT and CAW-5014.

**3. Materials:** The material used in the composition of the Latitude and Latitude EV Systems are:

**The humeral stem** is made from Cobalt-Chromium alloy according to ASTM F-799 Standard Specification for Cobalt-28 Chromium-6 Molybdenum Alloy Forgings for Surgical Implants (UNS R31537, R31538, R31539) (Latitude EV) or Cobalt-Chromium Molybdenum Casting Alloy according to ISO 5832-7 Implants for Surgery-Metallic Materials- Part 7 (Latitude);or Cobalt-chromium-molybdenum casting alloy according to ISO 5832-4 Implants for Surgery - Metallic Materials-Part 4 (Latitude). Only the Latitude EV Humeral stem is **coated** according to ASTM F1580 Standard Specification for Titanium and Titanium-6Aluminum-5Vanadium Alloy Powders for Coatings of Surgical Implants

**The Ulnar Stem** is made of Cobalt-28 Chromium-6 Molybdenum Alloy Forgings for Surgical Implants according to ASTM F-799 Standard Specification (UNS R31537, R31538, R31539) (Latitude EV) or Cobalt-Chromium Molybdenum Casting Alloy according to ISO 5832-7 Implants for Surgery-Metallic Materials- Part 7 (Latitude) or Wrought Cobalt-chromium-molybdenum alloy according to ISO 5832-12 Implants for Surgery–Metallic Materials–Part 12 (Latitude). Only the Latitude EV Ulnar stem is **coated** according to ASTM F1580 Standard Specification for Titanium and Titanium-6Aluminum-5Vanadium Alloy Powders for Coatings of Surgical Implants

**The humeral spool** is made of chromium-cobalt alloy (CrCo) according to ISO standard 5832-7 or ISO 5832-12 and of polyetheretherketone (PEEK-OPTIMA polymer) according to standard ASTM F 2026.

**The humeral screw** is made of stainless steel according to ISO 5832-9.

**The Ulnar radial components** are made of chromium-cobalt alloy (CrCo) according to ISO standard 5832-7 or ISO standard 5832-12 or ISO standard 5832-4 and ultra high molecular weight polyethylene (UHMWPE) according to ISO standard 5834-2.

**4. Indications:** The Tornier Elbow Prosthesis is intended for total elbow arthroplasty. Prosthetic replacement with this device may be indicated to relieve severe pain or significant disability following the effects of primary or secondary osteoarthritis and rheumatoid arthritis; correction of functional deformities; revision procedures where other treatments or devices have failed; treatment of fractures that are unmanageable using other techniques.

When the anatomic spool is used solely with the humeral stem, the Latitude elbow prosthesis is a distal hemi-elbow prosthesis indicated for the same indications that are previously described.\*\*

The prosthesis is for **cemented use** only.

**\*\*Note: The Anatomic humeral spool is not available for sale within the United States.**

**5. Known contra-indications to date:** Systemic infection is an absolute contraindication. Every effort should be made to rule out the possibility of preoperative sepsis in a patient who has one or more of the following abnormalities: fever and/or local inflammation; rapid joint destruction or bone resorption apparent on roentgenograms, and elevation of sedimentation rate unexplained by other disease, elevation of WBC count. Distant foci of infection from genitourinary, pulmonary, skin and other sites, dental focus infection which may cause hematogenous spread to the implant site. Skeletally immature patients. Cases where there is inadequate neuromuscular status, poor bone stock, or poor skin coverage around the elbow joint that would make the procedure unjustifiable. Neuromuscular or psychiatric disorders which might jeopardize fixation and postoperative care. Known allergy to one of the materials. Pregnancy.

**6. Side-effects and possible complications:** Dislocation, sepsis, loosening, instability, bone fracture, ulnar nerve palsy, neuropathy, triceps rupture, wound problems, components disassembly or breakage, and metal sensitivity in patients following total joint replacement have been rarely reported.

**7. Warnings and Cautions:**

The Latitude and Latitude EV Total Elbow Prosthesis has not been evaluated for safety and compatibility in the MRI environment. The system has not been tested for heating or migration in the MR environment.

Although published studies indicate that the MRI procedure has minimal effects on the material in this product, patients should be informed that there are several different manufacturers and generations of MRI equipment, and Tornier Inc. cannot make claims regarding the compatibility of this product with any specific MRI unit. It is recommended that patients and physicians consult with the MRI equipment manufacturers to discuss the compatibility of this product with the MRI equipment before undergoing any MRI procedure.

• This product is provided sterile in an undamaged package and is for ONE TIME USE ONLY. If either the implant or the package appears damaged, expiration date has been exceeded or if the sterility is questioned for any reason, the implant should not be used. DO NOT clean, re-sterilize or reuse as this may damage or compromise performance of the devices and may expose patient to risk of transmitting infectious diseases.

• **The humeral screw design for Latitude and the humeral screw design for Latitude EV are NOT interchangeable.**

• **In revisions, the screw selected should match the humeral stem version: Latitude Humeral Stems should be used only with Latitude Humeral Screws and Latitude EV Humeral Stems should only be used with Latitude EV Humeral Screws.**

**Preoperatively:** The surgeon must be fully conversant with all the aspects of either the Latitude or Latitude EV surgical techniques and know the indications and contra-indications of these types of implants. The surgeon must be acquainted before the operation with the specific operative technique as said above and also of the products, which is available from the manufacturer. As part of the pre-operative examination, the surgeon must check that no biological, biomechanical or other factors are present that will affect the correct conduct of the operation and the postoperative period. Use provided x-ray templates from manufacturer to assist examination. An appropriate range of sizes must be available at the time of the operation. Check that the implant has not suffered from any deterioration.

**Intraoperatively:** The correct selection of the type and size of the implant appropriate to the patient and the positioning of the implant are extremely important. The functional surfaces of the implants must not suffer from any damage, shock, abrasion or other deterioration.

In case of change of the humeral spool, never re-use a humeral spool and/or screw that were already assembled. Also be sure that the appropriate compatible humeral screw is used for the appropriate Latitude or Latitude EV system.

In case of revision of the ulnar plating bushing, the ulnar cap **MUST** be assembled to the ulnar stem in order to secure the bushing.



**Postoperatively:** Patients must be informed by the surgeon of the precautions to be taken during daily life activities to guarantee the maximum implant service life. It is recommended that a regular postoperative follow-up be undertaken to detect early signs of wear, loosening of the prosthesis, etc., and to consider the action to be taken. Normal wear of the implant with respect to the state of knowledge at the time of its design cannot in any way be considered to constitute a dysfunction or a deterioration in the characteristics of the implant. A suitable rehabilitation program must be designed and implemented. Failure of the patient to comply with postoperative care instructions and rehabilitation may compromise surgical outcomes.

**8. Storage and handling:** Implants must be stored in their original sealed packaging. Storage place must be away from humid environments. Implants must not be exposed to direct sunlight, ionizing radiation, extreme temperatures nor particular contamination. Implants must be handled with care to preserve integrity of their packaging.

**9. Packaging and sterilization:** The implants are supplied sterile (gamma radiation). The expiration date for sterilization must be checked prior to surgery. Only products implanted before the expiration date may be considered sterile. The packaging and labeling must be checked for integrity. Reject any implant if the packaging shows any flaw. Every precaution must be taken to ensure sterility when opening the packaging of the implant and when inserting it.

**LIMITED WARRANTY**

Tornier, Inc. warrants that this product meets the manufacturer’s specifications and is free from manufacturing defects at the time of delivery. This warranty specifically excludes defects resulting from misuse, abuse, or improper handling of the product subsequent to receipt by the purchaser. Tornier, Inc. does not warrant the outcome of the surgical procedure.

<b>Symbols Key:</b>		
	CATALOG NUMBER	 DO NOT REUSE
	BATCH CODE	 DO NOT USE IF PACKAGE IS DAMAGED
	SERIAL NUMBER	 USE BY
	AUTHORIZED REPRESENTATIVE OF THE EUROPEAN COMMUNITY	 CAUTION: FEDERAL LAW (USA) RESTRICTS THIS DEVICE TO SALE BY OR ORDER OF A PHYSICIAN
	STERILIZED USING IRRADIATION	
	CAUTION	 MANUFACTURER
	CONSULT INSTRUCTIONS FOR USE	

**CONTACT INFORMATION**

**Manufacturer:**

Tornier, Inc.  
10801 Nesbitt Avenue South  
Bloomington, MN 55437  
USA  
Tel: (952) 921-7100  
Fax: (952) 236-4007

**European Authorized Representative and Distributor outside of USA and Canada:**

 Tornier S.A.S  
161 Rue Lavoisier  
38330 Montbonnot Saint Martin  
France  
+33(0)4 76 61 35 00  
Fax: +33(0)4 76 61 35 33

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©Copyright 2015 Tornier, Inc.