ENGLISH (EN)

# TORNIER

## Force Fiber®

POLYETHYLENE SUTURE

Non-absorbable Surgical Suture

Meets USP except for oversized diameter.

White Braided - Ultra High Molecular Weight Polyethylene

White/Blue Co-braid - Ultra High Molecular Weight Polyethylene

White/Black Co-braid - Ultra High Molecular Weight Polyethylene

#### DESCRIPTION

Force Fiber® White Braided Polyethylene Surgical Suture is a non-absorbable, sterile surgical suture composed of undyed (white) ultra high molecular weight polyethylene.

Force Fiber® White/Blue Co-braid is a non-absorbable, sterile surgical suture composed of undyed (white) ultra high molecular weight polyethylene braided with one of two strands of blue polypropylene to add color. The blue colorant is [phthalocyaninato(2-)] copper with a concentration not to exceed 0.5% by weight.

Force Fiber® White/Black Co-braid is a non-absorbable, sterile surgical suture composed of undyed (white) ultra high molecular weight polyethylene braided with one strand of black polyamide to add color. The black colorant is a lewood extract (hematein dye) with a concentration not to exceed 1.0% by weight.

Force Fiber® sutures exceed USP specifications for diameter.

Suture Size	USP Ave. Diameter Specification (mm) <861>	Maximum Oversize Average Diameter (mm)	Maximum Oversize Average Diameter from USP (mm)
2	0.500 - 0.599	0.630	0.031
3&4	0.600 - 0.699	0.754	0.055
5	0.700 - 0.799	0.861	0.062

#### INDICATIONS

Force Fiber® Non-absorbable Surgical Sutures are indicated for use in approximation and/or ligation of soft tissues, including use of allograft tissue for orthopedic surgeries.

#### CONTRAINDICATIONS

None known (US).

Force Fiber® Non-absorbable Surgical Sutures are contraindicated for cardiovascular uses (EU).

#### WARNINGS

■ Do not use if sealed package is tampered with or damaged.
■ Discard sutures that have been opened even if unused.
■ Suture should NOT be Re-sterilized.

Force Fiber® Non-absorbable Surgical Suture may elicit a minimal acute inflammatory reaction in tissues, which is followed by gradual encapsulation of the suture by fibrous connective tissue.

Users should be familiar with surgical procedures and techniques involving non-absorbable sutures before employing Polyethylene Surgical Suture for wound closure, as the risk of wound dehiscence may vary with the site of application and the suture material used.

Acceptable surgical practice must be followed with respect to drainage and closure of infected or contaminated wounds.

Re-processing of products intended for single use only may result in degraded performance or a loss of functionality. Re-use of single use only products may result in exposure to viral, bacterial, fungal, or prionic pathogens. Validated cleaning and sterilization methods and instructions for reprocessing to original specifications are not available for these products. This product is not designed to be cleaned, disinfected, or sterilized.

#### PRECAUTIONS

Avoid excessive handling i.e. crushing or crimping resulting from use of surgical instruments such as forceps and needle holders.

Adequate knot security requires the accepted surgical technique of flat, square ties, with additional throws as warranted by surgical circumstance and the experience of the surgeon.

Strong familiarity with surgical procedures and techniques for non-absorbable sutures is required before use.

Discard used needles in appropriate container ("sharps"). This device should be handled and disposed of in accordance with all applicable regulations including, without limitation, those pertaining to human health and safety and the environment.

#### ADVERSE REACTIONS

Adverse effects associated with the use of this device include: wound dehiscence, calcuI formation in urinary and biliary tracts when prolonged contact with salt solutions such as urine and bile occurs, infected wounds, minimal acute inflammatory tissue reaction, and transitory local irritation.
As with any foreign body, prolonged contact of this or any other suture with salt solutions, such as those found in urinary or biliary tracts, may result in calculus formation.


#### HOW SUPPLIED

Force Fiber® Surgical Suture is provided sterile and braided;
■ Undyed white polyethylene.
■ Co-braid of undyed white polyethylene with blue monofilament polypropylene.
■ Co-braid of undyed white polyethylene with black monofilament nylon (polyamide).

The suture is available in USP sizes 2, 3&4, and 5 (metric sizes 5, 6, and 7). The suture is provided in a variety of lengths with and without needles. Finished suture may be packaged in cartons as single packs, multi-packs, or procedure packs.

	Catalog Number
	Batch Code
	Authorized Representative in The European Community
	Manufacturer
	Do not use if package has been opened or damaged
	Single use only
	Do not Re-sterilize
	Sterilized using ethylene oxide
	Use by
	Consult instructions for use
	Caution

Force Fiber® is a registered trademark of Teleflex Incorporated.
TORNIER® is a registered trademark of Tornier, Inc.

	Teleflex Medical 2917 Weck Drive Research Triangle Park, NC 27709 USA International: (919) 544-8000    USA: (866) 246-6990
	Teleflex Medical IDA Business and Technology Park Dublin Road, Athlone, Co. Westmeath, Ireland
	Distributed by: Tornier, Inc. 10801 Nesbitt Avenue South Bloomington, MN 55437 USA Tel: +1 952-921-7100 Fax: +1 952-236-4007
	Tornier S.A.S. 161 Rue Lavoisier 38330 Montbonnot Saint Martin, France Tel: +33 (0)4 76 61 35 00 Fax: +33 (0)4 76 61 35 33
	0120 www.tornier-us.com
	©Copyright 2013 Teleflex Medical

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TORNIER, Inc. 10801 Nesbitt Avenue South Bloomington, MN 55437 USA Tel: +1 952-921-7100 Fax: +1 952-236-4007