## BIOFIBER™ CM SCAFFOLD Surgical Mesh

IFU-5171 Rev. C

Please read this information carefully before using the BioFiber™ CM Scaffold

### DESCRIPTION

BioFiber™ CM Scaffold is a bi-layer, synthetic absorbable reinforced woven fabric made from poly-4-hydroxybutyrate fibers coated with bovine-derived collagen. The construction permits the mesh to be cut into any desired shape or size with minimal unraveling. The device is supplied sterile, for single-patient use for the reinforcement of soft tissue where weakness exists.

### INDICATIONS

BioFiber™ CM Scaffold is intended for reinforcement of soft tissues, in conjunction with sutures and/or suture anchors during tendon repair surgery; including the reinforcement of rotator cuff, patellar, Achilles, biceps and quadriceps tendons.

BioFiber™ CM Scaffold is not intended for use as a replacement for normal body structures or to provide the full mechanical strength to support the rotator cuff, patellar, Achilles, biceps, or quadriceps tendons. Sutures, used to repair the tear, and sutures or bone anchors, used to attach the tissue to the bone, provide mechanical strength for the tendon repair.

### CONTRAINDICATIONS

Because BioFiber™ CM Scaffold is absorbable; it should not be used where permanent wound or organ support is required.

The device is contraindicated for use in patients with:

- Active or latent infection,
- Decreased vascularity,
- Pathologic soft tissue conditions

The device is also contraindicated in uses that require rolling, folding, or layering, and which may create a space impermeable to fluid, cells, and blood vessels. Such uses may result in excessive inflammation, drainage, extrusion or infection.

### WARNINGS

Use of this product in applications other than those intended for implantation to reinforce soft tissue where weakness exists has the potential for serious complications. The patient should be made aware of the potential complications as listed.

This product should not be used in individuals who are sensitive or allergic to bovine collagen.

Device manufacture involves exposure to tetracycline hydrochloride and kanamycin sulfate. The safety and product use for patients with hypersensitivities to these antibiotics is unknown.

The safety and effectiveness of BioFiber™ CM Scaffold in neural tissue and in cardiovascular tissue has not been established.

### PRECAUTIONS

- Do not use this product without reading and understanding these instructions.
- Do not resterilize. Discard all open and unused portions of BioFiber™ CM Scaffold.
- The device is sterile if the package is unopened and undamaged. Do not use if the package seal is broken.
- Discard the device if mishandling has caused possible damage or contamination.
- Aseptic technique must be adhered to throughout the procedure.
- Discard BioFiber™ CM Scaffold if past Use-by-date of the product [as indicated as year (4 digits) and month (2 digits)]. The product is not to be used after the last day of the month indicated.
- Rectangular patches smaller than 30mm x 30mm and circular patches smaller than 20mm diameter should not be trimmed.
- Trimming patches will decrease suture retention strength. If trimming of patches is necessary, it is recommended to place mattress stitches no less than 5 mm from cut edge. Anticipate a 20% reduction in suture retention strength in doing so.
- DO NOT clean, resterilize, or reuse as this may damage or compromise performance of devices and may expose patient to risk of transmitting infectious disease.

### ADVERSE EFFECTS

The following potential complications are possible with the use of surgical mesh material. If any of these conditions occur, the device may be removed at the surgeon's discretion.

- Infection
- Wound dehiscence.
- Acute or chronic inflammation.
- Allergic reactions to device materials
- Tissue erosion
- Device breakage or migration.

### ACTIONS

BioFiber™ CM Scaffold degrades through a process of hydrolysis and hydrolytic enzymatic digestion. Implantation studies indicate that the absorption of the mesh material is minimal until about 26 weeks post implantation and essentially complete in 52 weeks.

### DIRECTIONS FOR USE

1. Prepare the implant site using standard surgical techniques.
2. Using aseptic technique, trim the mesh to fit the implant site, providing a small allowance for overlap. See PRECAUTIONS section for restrictions on trimming patches and potential impact on strength.
3. Transfer the scaffold to the graft site and suture into place, avoiding excess tension.
4. Secure the scaffold by placing sutures 5 – 10mm apart and 5mm from the edges of the mesh.
5. Complete the standard surgical procedure.
6. Discard any unused portions of the scaffold.

### STERILITY

BioFiber™ CM Scaffold is sterilized by ethylene oxide gas (EO). Product is sterile unless package has been opened or damaged. DO NOT RESTERILIZE. Discard opened or unused containers.

### HOW SUPPLIED

The BioFiber™ CM Scaffold is available in a variety of dimensions and may be cut to size as required for the specific use.

### STORAGE AND HANDLING

Store in a cool, dry place and in a manner that protects the integrity of the package and sterile barrier.

### PACKAGING AND LABELING

- Tornier BioFiber™ CM Scaffold should be accepted only if the factory packaging and labeling arrive intact.
- Contact Customer Service if package has been opened or altered.

**CAUTION:** Federal law (USA) restricts this device to sale by or on the order of a physician.

### SYMBOL KEY

	CATALOG NUMBER		NOT MADE WITH NATURAL RUBBER LATEX
	BATCH CODE		DO NOT USE IF PACKAGE IS DAMAGED
	STERILISED USING ETHYLENE OXIDE		TEMPERATURE LIMITATION
	CAUTION, CONSULT ACCOMPANYING DOCUMENTS		USE BY
	CONSULT INSTRUCTIONS FOR USE		CAUTION: FEDERAL LAW (USA) RESTRICTS THIS DEVICE TO SALE BY OR ORDER OF A PHYSICIAN
	DO NOT RESTERILISE		MANUFACTURER
	DO NOT REUSE		

### CUSTOMER SUPPORT

For product complaints and potential adverse events, please contact your local Sales Representative or Customer Support.

**Tornier, Inc.**  
10801 Nesbitt Avenue South  
Bloomington, Minnesota  
55437 USA  
Tel: +1 952 921 7100  
Fax: +1 952 236 4007

© Copyright 2015 Tornier, Inc.

Tornier® is a registered trademark of Tornier, Inc. BioFiber™ and BioFiber™ CM are also Tornier Trademarks. These products are covered by one or more U.S. and international patents pending. All rights reserved. Printed in the USA.

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