FRANCAIS

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trapeziometacarpal joint without osteoarthritis

- preserve the trapezium, thereby providing an alternative solution, less invasive than trapezectomy
- consequently preserve the length of the thumb column, and therefore its kinematics and thus painless mobility of it.

4. INDICATIONS

The CMI implant can be used:

- For painful rhizarthrosis that is resistant to conservative treatment, for which the ligament system of the joint can ensure its stability with a ligamentoplasty
- Revision (this indication is not cleared in United States):
 - For partial trapezectomy or partial trapezectomy with or without ligamentoplasty
 - For trapeziometacarpal interposition or prosthesis
 - For trapeziometacarpal ligamentoplasty

- Revison (this indication is not cleared in United States):

- For partial trapezectomy or partial trapezectomy with or without ligamentoplasty
- For trapeziometacarpal interposition or prosthesis
- For trapeziometacarpal ligamentoplasty

Regardless of the radiological stage, using this implant requires a stabilization plasty at the expense of the extrinsic tendons of the 1st radius. Whether primary or post-traumatic, isolated trapeziometacarpal osteoarthritis can be treated with a CMI implant. In the event of peri-trapezial impairment, another type of implant or combination of treatments must be discussed.

5. KNOWN CONTRAINDICATIONS TO DATE

The use of the CMI implant is contraindicated in the following cases (non-exhaustive list):

- Sequelae of infections of the wrist,
- Major stiffness with muscle and tendon retraction or major instability of the 1st radius,
- Acute or chronic infectious diseases, regardless of the etiology and location,
- Neuromuscular or psychiatric diseases, which may cause the surgery to fail
- Inflammatory context or presence of components or debris from a previous prosthesis
- Bone immaturity - Growing patient with active epiphysis.
- Unsuitable or insufficient bone support not enabling the correct fitting of the implant;
- Inadequate condition of the skin, bones or neurovascular system
- Patients with high functional demand, especially: athlete, manual laborer,
- Known or suspected allergy to the material.
- Pregnant women

6. SIDE EFFECTS AND POSSIBLE COMPLICATIONS

- Migration or dislocation of the implant
- Subsidence of the implant into the bone.
- Functional changes of the anatomical components of the other

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surrounding bones

- Complications associated with wear of the component
- Rupture of the component subjected to abnormal load or activity
- Dimensional modification of the limb operated on
- Reduction in the degree of movement.
- Hematoma in the operated area
- Infection in the operated area
- Neuropathies
- Cardiovascular problem
- Tissue reaction
- Pain for several days

Certain complications and side effects may result from ignorance of the attached precautions for use.

7. PRECAUTIONS FOR USE

- Before deciding on proceeding with the implant, evaluate all other possible, more conservative, treatments to ensure that all such treatments will not provide results as good as those of the CMI.
- Never reuse an implant, even in seemingly good condition, in order to avoid any risk of cross-contamination or change in performance.
- Never re-sterilize an implant delivered sterile.

The following situations hinder the success of an interposition implant:

- Sports activities or high activity levels,
- People likely to fall,
- Alcoholism or drug use,
- Inability of the patient to follow the surgeon's recommendations and the rehabilitation program.

Preoperatively:

The surgeon must be fully conversant with all aspects of the surgical technique and know the indications and contra-indications for this type of implant. The surgeon must be acquainted, prior to surgery, with the specific operative technique of the product, available from the manufacturer. As part of the pre-operative examination, the surgeon must check that no biological, biomechanical or other factors are present, that will affect the correct conduct of the operation and the postoperative period. He must also check that the quality of the bone is satisfactory enough to support the implantation. An appropriate range of sizes must be available at the time of the operation.

Intraoperatively:

Correct selection of appropriate implant type and size for a patient, as well as implant position, is extremely important. The use of trial pieces allows for the proper size selection of the implants. Frequent radioscopic checks allow the position of the prosthesis to be checked.

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Handling of the implant with metal clamps is strictly not recommended.

The implants must not be used if their functional surfaces have been damaged or have undergone shock, abrasion, or other deterioration.

Post-operatively:

A suitable rehabilitation program must be designed and implemented specific to the patient. The surgeon must inform patients about:

- Precautions to take in daily life to guarantee maximum implant survival,
 - The fact that their weight and level of activity can affect the life span of the prosthesis,
- The implant has not been evaluated for safety and compatibility in the MR environment. The implant has not been tested for heating or migration in the MR environment.

It is recommended that a regular postoperative follow-up is undertaken to detect early signs of wear and to consider the action to be taken. Normal wear of the implant in respect of the state of knowledge at the time of its design cannot in any way be considered to constitute a dysfunction or deterioration in the characteristics of the implant.

The stability of the implant is ensured by its shape, but also by the preparation by the surgeon during fitting, and especially by the healthy condition of the bones and all the soft parts surrounding the implant (ligaments, tendons, muscles, etc.). Migration/dislocation of the implant can in no way be considered a malfunction or alteration of these characteristics.

8. STORAGE AND HANDLING

Implants must be stored in their original sealed packaging. Implants must be stored away from moisture. Implants must not be exposed to direct sunlight, ionizing radiation, extreme temperatures or particulate contamination. Implants must be handled carefully to preserve the integrity of their packaging.

9. PACKAGING AND STERILIZATION

Implants are sterilized by gamma radiation. The expiration date of sterilization and package integrity must be verified.

An implant whose packaging is open or damaged or whose expiration date has passed must not be used.

All aseptic precautions must be taken when opening the sterile packaging of the implant and during implantation. Ancillary instruments are supplied non-sterile. For handling and sterilization of non-sterile ancillary instruments, refer to the instruments instructions for use.

The templates are supplied non-sterile and should not be sterilized. For any other information regarding the ancillary instruments, refer to the instructions provided for this purpose.

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Hersteller: TORNIER SAS
Geschäftsanschrift in den USA: TORNIER, Inc