## FRANCAIS

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progression in the diaphysis and therefore facilitate the cement pressurization during the introduction of the implant

**4. Indications for use:**

The Tornier Cement Restrictor is a diaphyseal plug designed to occlude the medullary cavity before the introduction of acrylic cement.

**5. Instructions for use:**

Before cementing the stem, the Tornier Cement Restrictor is inserted into the bone shaft by using the inserter. Three sizes are available:

- the small size (dia. 13 mm) is indicated for diaphyseal diameters from 5 mm to 7 mm after boring,
- the medium size (dia. 24 mm) is indicated for diaphyseal diameters from 7 mm to 15 mm after boring,
- the big size (dia. 38 mm) is indicated for diaphyseal diameters from 12 mm to 20 mm after boring.
The insertion depth is calculated by adding 10 mm to the stem length. Once the plug is positioned into the shaft, the inserter is removed and the cement is injected.

**6. Known contra-indications to date:**

- Acute or chronic infectious diseases of any aetiology and localization,
- unsuitable or insufficient bone support,
- bone immaturity,
- known allergy to the material.

**7. Side-effects and possible complications:**

- Infection.
- Break of wing.
- Some complications and side-effects may stem from a lack of awareness of the precautions for use mentioned below.

**8. Warnings and cautions:**

- **Never re-use an implant, even if it seems to be in perfect condition, to prevent any risks of cross-contamination or a risk of reduced performances.**

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implant has not been tested for heating or migration in the MR environment.

There is no proven incompatibility of the prosthesis to the MR environment. However, as a precaution, it is recommended not to expose a patient wearing the prosthesis to the MR environment.

It is recommended that a regular postoperative follow-up is undertaken to detect early signs of wear, loosening of the implant, etc., and to consider the action to be taken. A suitable rehabilitation program must

**9. Storage and handling:**

Implants must be stored in their original sealed packaging. Storage place must be away from humidity. Implants must not be exposed to direct sunlight, ionising radiation, extreme temperatures nor particular contamination. Implants must be handled with care to preserve integrity of their packaging.

**10. Packaging and sterilization:**

The implants are supplied sterile (gamma radiation). The expiration date for sterilization and integrity of the packaging must be checked.

An implant whose packaging is open or damaged or whose expiration date has passed must not be used.

Every precaution must be taken to ensure sterility when opening the packaging of the implant and during implantation.

Ancillary instruments may be supplied sterile.




For handling and sterilization of non-sterile ancillary instruments, refer to the ancillary instruments instructions. The templates are supplied non-sterile and should not be sterilized.

**11. Implant Retrieval and handling:**

In case of retrieval of the implant from the patient; the retrieved implant should be handled according to appropriate and validated hospital procedures.

For any other information regarding the ancillary instruments, refer to the instructions provided for this purpose.

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


**Interpretation of terms and symbols**

|   |   |
|---|---|
|    | Do not reuse                            |
|    | Expiration date                         |
|    | Sterilized using irradiation            |
|    | Batch code (sterilization)              |
|   | Serial number (manufacturing)           |
|  | Catalogue number                        |
|  | Caution, consult accompanying documents |
|  | Consult instructions for use            |
|  | Do not use if package is damaged        |
|  | Do not re-sterilize                     |
|  | Not made with natural rubber latex      |

The French text is the reference text.


CE marking is only valid if it is also mentioned on the external package labeling

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**Hersteller: TORNIER SAS**
**TORNIER, Inc**
**Fabrikant: TORNIER SAS**