



OSTEOSET[®] Resorbable Bead Kits and Mini-Bead Kits Mixing Instructions

135040-3

The following languages are included in this packet:

English (en)

Deutsch (de)

Nederlands (nl)

Français (fr)

Español (es)

Italiano (it)

Português (pt)

中文-Chinese (sch)

Türkçe (tk)

For additional languages, visit our website www.wmt.com


Then click on the **Prescribing Information** option.

For additional information and translations please contact the manufacturer or local distributor.



CE 0086*

Wright Medical Technology, Inc.
P.O. Box 100
Arlington, TN 38002
U.S.A.



Wright Medical EMEA
Krijgsman 11
1186 DM Amstelveen
The Netherlands

*** The CE-Marking of Conformity is applied per catalog number and appears on the outer label, if applicable.**

Rx ONLY

November 2009
Printed in U.S.A.

MIXING INSTRUCTIONS FOR THE OSTEOSSET[®] RESORBABLE BEAD KIT (135040-3)

1. Mix according to directions listed in Steps 1-5.

Note: Do not add any additional substances to the paste. Using alternative mixing solutions and/or adding other substances to the mixture may alter the setting time significantly. Some substances, such as bone marrow and blood, will prevent the paste from setting altogether.



Step 1

- Add powder to mixing bowl
- Add all of diluent
- Allow to sit
 - **Standard kit: 1 minute**
 - **Fast-cure kit: 30 seconds**

Step 2

- Mix thoroughly for 30 – 45 seconds
- DO NOT over mix
- When the consistency of the mixture is paste-like (i.e. sticks to the spatula), it is ready to apply to mold




Step 3

- With spatula apply an even coat of calcium sulfate paste to the Wright Bead Template bottom half (shown), ensuring complete filling of each bead cavity. Each OSTEOSSET[®] Resorbable Bead Kit contains enough material to fill two bead templates.



Step 4

- Apply mold top and firmly press together (excess material will extrude from vents on the top of the mold). Allow to sit:
Standard Kit: approximately 15-20 minutes
Fast Cure Kit: approximately 3-5 minutes
- Carefully separate the mold



Step 5

- Flex the mold to facilitate bead removal

2. Gently pack the molded implants into the treatment site. Avoid overfilling the bone void or compressing the treatment site.
3. Remove excess material from the treatment site. Cover the opening of the bone void with cortical bone window and/or healthy periosteum or other soft tissues. Close the site using standard closure techniques.
4. Discard any unused OSTEOSSET[®] Resorbable Bead Kit material and the mixing tools.

MIXING INSTRUCTIONS FOR OSTEOSSET[®] MINI RESORBABLE BEAD KIT

1. Mix according to directions listed in Steps 1-5.

Note: Do not add any additional substances to the paste. Using alternative mixing solutions and/or adding other substances to the mixture may alter the setting time significantly. Some substances, such as bone marrow and blood, will prevent the paste from setting altogether.

Step 1

- Add powder to mixing bowl
- Add all of diluent
- Allow to sit
 - **Standard kit: 1 minute**
 - **Fast-cure kit: 30 seconds**




Step 2

- Mix thoroughly for 30 – 45 seconds
- DO NOT over mix
- When the consistency of the mixture is paste-like (i.e. sticks to the spatula), it is ready to apply to mold



Step 3


- With spatula apply an even coat of OSTEOSSET[®] paste to only one half of the supplied mold
- Ensure complete filling of each bead cavity.
NOTE: Each kit contains enough material to fill only one half of the supplied mold, which totals approximately 5cc final volume.
- After completely filling the mold half, tap on a flat surface to remove air bubbles.



Step 4


Allow to sit:

- **Standard Kit: approximately 15-20 minutes**
- **Fast Cure Kit: approximately 3-5 minutes**



Step 5

- Flex the mold to facilitate bead removal



2. Gently pack the molded implants into the treatment site. Avoid overfilling the bone void or compressing the treatment site.
3. Remove excess material from the treatment site. Cover the opening of the bone void with cortical bone window and/or healthy periosteum or other soft tissues. Close the site using standard closure techniques.
4. Discard any unused OSTEOSSET[®] Mini Resorbable Bead Kit material and the mixing tools.