
Distributed By:  Wright Medical Technology, Inc. 901.867.9971   800.238.7188   www.wmt.com	Manufactured By:  11621 Research Circle, Alachua FL USA 32615 386.418.8888   800.624.7238   www.rtx.com
---	--

#### CANCELLO-PURE® COTTON, EVANS, 10 X 50MM WEDGE

	Read this entire package insert carefully prior to use.		Single patient use only, on a single occasion.		Restricted to sale by or on the order of a licensed healthcare practitioner (e.g., physician).
---	---	---	--	---	--

CANCELLO-PURE® implants are bovine tissue procured and processed according to U.S. Department of Agriculture and U.S. Food & Drug Administration regulations. This implant has been processed to meet high levels of quality and safety, and has passed quality inspection criteria. The implant is provided sterile. Supplemental fixation is required for use in osseous defects that are intrinsic to the stability of the bony structure.

##### INDICATIONS FOR USE

These implants are indicated for bony voids or gaps that are not intrinsic to the stability of the bony structure. They are indicated to be placed into bony voids or gaps of the skeletal system (e.g., extremities, spine, ilium and/or pelvis). These defects may be surgically created osseous defects or osseous defects created from traumatic injury to the bone. The implant provides a void filler that remodels into the recipient’s skeletal system.

##### CONTRAINDICATIONS

There are no contraindications for this implant.

##### WARNINGS AND PRECAUTIONS

The same medical/ surgical conditions that may complicate any surgical procedure may occur during or following the transplantation of this implant. Although *in vitro* and *in vivo* testing showed no significant antigenic response to this implant, there is a potential for immune response following implantation.

##### STORAGE

Store at temperature range specified on the product label.

##### PROCESSING

	BioCleanse® is a validated sterilization process that inactivates or removes potential pathogens through a complex, proprietary combination of chemical treatments and mechanical processes.
	Low temperature, low dose gamma irradiation is applied terminally to the product to achieve a sterility assurance level (SAL) of 10 <sup>-6</sup> .

Trace amounts of manufacturing residuals may remain after processing; acetic acid, detergents, hydrogen peroxide and isopropyl alcohol.

##### GENERAL INSTRUCTIONS

- Use on a single occasion for a single patient only. Once the package is opened, the implant must be used for the current procedure or discarded.
- The outermost packaging is non-sterile and is used to protect the implant during shipping and storage.
- The product’s sterile barrier packaging is comprised of two sealed pouches. To prevent contamination of the implant, use sterile technique for preparation and implantation.
- Remove the triple barrier packaged product and package insert from the box.
- Inspect the implant, packaging and label materials carefully:
  - Do not use past expiration date.
  - Do not use if the implant or packaging is damaged.
  - Do not use if there are discrepancies in label information.
- Additional product should be available in case of an unexpected need during the procedure.
- Do not re-sterilize the implant.
- Discard any unused implant in accordance with applicable facility protocol.
- Promptly report all product defects and patient adverse events to Wright Medical Technology, Inc. (See Customer Returns & Complaints section).

##### IMPLANT PREPARATION

- Open the outer package and pass the sterile inner package to the sterile field.
- Open both layers of the inner package within sterile field.
- Place implant into sterile bowl and cover with sterile water or physiological saline. At the surgeon’s discretion, the following autologous fluids; blood, bone marrow aspirate and platelet rich plasma, may also be used to hydrate the implant.
- Hydrate for a minimum of 30 seconds before use.

##### WARRANTY STATEMENT

*This biologic graft, processed and packaged for surgical implantation, is unique and does not constitute a product under liability laws. No implied warranties of merchantability or fitness for a particular purpose are applicable. No implied warranties exist as to defects in biologics which cannot be detected, removed, or prevented by reasonable use of available scientific procedures or techniques. Furthermore, ALL WARRANTIES ARE DISCLAIMED, WHETHER EXPRESSED OR IMPLIED BY OPERATION OF LAW OR OTHERWISE INCLUDING ALL IMPLIED WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE.*

##### CUSTOMER RETURNS & COMPLAINTS


For further information, or to report a complaint or adverse reaction:

<b>For U.S.A. contact:</b> Wright Medical Technology, Inc. 5677 Airline Road Arlington, TN 38002, USA Tel: 901.867.9971   800. 238.7188	<b>For Europe contact:</b> Wright Medical EMEA Krijgsman 11 1186 DM Amstelveen The Netherlands Tel: +31 (0) 20.545.01.00	<b>EC REP</b> Emergo Europe Molenstraat 15 2513 BH The Hague The Netherlands Tel: (+31) 70.345.8570 Fax: (+31) 70.346.7299
---	---	---

CANCELLO-PURE is a registered trademark of Wright Medical Technology, Inc.

BioCleanse is a registered trademark of RTI Biologics, Inc.

Note: The CE-Marking of conformity is applied per catalog number and appears on the outer label, if applicable.

DEFINITION OF LABEL SYMBOLS			
	Caution, see instructions for use		Sterile by Gamma Irradiation
	LOT		Lot number
	REF		Single patient use (Do not reuse)
	Storage temperature limits		Authorized Representative in the European Community
			Expiration Date

For use by RTI and affiliates. This document contains product specific information that may not be reproduced without prior written approval.

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CANCELLO-PURE ist eine eingetragene Marke von Wright Medical Technology, Inc.
BioCleanse ist eine eingetragene Marke von RTI Biologics, Inc.