






















Symbols Glossary

Symbol	Symbol Title	Symbol Description	Standard Reference	Title and Designation Number of the Standard
	Manufacturer	Indicates the medical device manufacturer, as defined in EU Directives 90/385/EEC, 93/42/EEC and 98/79/EC.	5.1.1	ISO 15223-1:2016 Medical Devices – Symbols to be used with medical device labels, labelling and information to be supplied – Part 1: General Requirements
	Authorized representative in the European Community	Indicates the authorized representative in the European Community.	5.1.2	
	Date of manufacture	Indicates the date when the medical device was manufactured.	5.1.3	
	Use-by date	Indicates the date after which the medical device is not to be used.	5.1.4	
	Batch code	Indicates the manufacturer's batch code so that the batch or lot can be identified.	5.1.5	
	Catalogue Number	Indicates the manufacturer's catalogue number so that the medical device can be identified.	5.1.6	
	Sterile	Indicates a medical device that has been subjected to a sterilization process.	5.2.1	
	Sterilized using aseptic processing techniques	Indicates a medical device that has been manufactured using accepted aseptic techniques.	5.2.2	
	Sterilized using ethylene oxide	Indicates a medical device that has been sterilized using ethylene oxide.	5.2.3	
	Sterilized using irradiation	Indicates a medical device that has been sterilized using irradiation.	5.2.4	
	Do not re-sterilize	Indicates a medical device that is not to be re-sterilized.	5.2.6	
	Non-sterile	Indicates a medical device that has not been subjected to a sterilization process.	5.2.7	
	Do not use if package is damaged	Indicates a medical device that should not be used if the package has been damaged or opened.	5.2.8	
	Keep away from sunlight	Indicates a medical device that needs protection from light sources.	5.3.2	
	Keep dry	Indicates a medical device that needs to be protected from moisture.	5.3.4	
	Temperature limit	Indicates the temperature limits to which the medical device can be safely exposed.	5.3.6	
	Do not re-use	Indicates a medical device that is intended for one use, or for use on a single patient during a single procedure.	5.4.2	
	Consult instructions for use	Indicates the need for the user to consult the instructions for use.	5.4.3	
	Caution	Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions that cannot, for a variety of reasons, be presented on the medical device itself.	5.4.4	
	MR Conditional	An item with demonstrated safety in the magnetic resonance environment within defined conditions.	N/A	
	Prescription Only	Caution: U.S. Federal Law restricts this device to sale by or on the order of a physician	N/A	N/A

