

SIDEKICK® EZ FRAME™

External Fixator System



SIDEKICK® EZ FRAME™

External Fixator System

External Fixation, Simplified

The SIDEKICK® EZ FRAME™ was designed to address the need for a simplified approach to external fixation. The system allows for the management of diabetic Charcot reconstruction, triple arthrodesis, fusion of the midfoot and rearfoot, challenging comminuted trauma, and most foot pathology not requiring fixation above the ankle.



By utilizing a unique combination of total contact casting concepts for the tibia and thin wire fixation in the foot, tibial wires or pins are not needed. By requiring no pins or wires in the tibia, the SIDEKICK® EZ FRAME™ eliminates complications associated with tibial fixation elements and the time it takes to place these pins and wires.

Advanced Design Features

- Soft foot pad
- Customizable liner
- Adjustable rocker-bottom plate
- Detachable rocker-bottom sole

Key System Benefits

- Simplify frame application
- Eliminate tibial pins and wires
- Off-load the plantar aspect of the foot as needed
- Designed for enhanced patient comfort



Wright Medical Technology, Inc.
1023 Cherry Road
Memphis, TN 38117
800 238 7117
901 867 9971
www.wmt.com

Wright Medical EMEA
Atlas Arena, Australia Building
Hoogoorddreef 7
1101 BA Amsterdam
the Netherlands
011 31 20 565 9060