


PRO-TOE[®] C2

Hammertoe Fixation System




PRO-TOE® C2

Hammertoe Fixation System



System Features

- Addresses hammertoe correction and MTPJ instability with cannulated implants
- Unique insertion method provides compression across PIPJ
- Stainless steel implant offering includes 4 diameters in 0° angulation



Implant Offering

Proximal Diameter	Proximal Length	Distal Diameter	Distal Length	Style	Accepted K-wire
2.0mm	13mm	4.0mm	6.0mm	Solid	–
2.4mm	13mm	4.0mm	6.0mm	Solid	–
3.0mm	13mm and 16mm	4.0mm	6.0mm	Cannulated	1.1mm
3.5mm	13mm and 16mm	4.0mm	6.0mm	Cannulated	1.1mm and 1.6mm



1023 Cherry Road
 Memphis, TN 38117
 800 238 7117
 901 867 9971
www.wright.com

62 Quai Charles de Gaulle
 69006 Lyon
 France
 +33 (0)4 72 84 10 30
www.tornier.com

Unit 1, Campus Five
 Letchworth Garden City
 Hertfordshire SG6 2JF
 United Kingdom
 +011 44 (0)845 833 4435